

MEDICAL MIRoTEC is a therapeutic product designed to assist in the treatment of various musculoskeletal conditions and injuries.

It utilizes "space blanket" technology, consisting of a micro-perforated, metallized reflective foil laminated to a supportive polyester substrate.

When applied to the body, MIRoTEC reflects infrared radiation and enhances blood circulation while allowing moisture to pass through.

The increased blood circulation promoted by MIRoTEC helps facilitate healing in injuries, reduces inflammation in soft tissue and joint disturbances, and alleviates stiffness in the back and muscle strains.

This improved circulation can be particularly beneficial for individuals with arthritis, enhancing mobility and providing pain relief.


MIRoTEC can also be used proactively to prevent injuries and aid in recovery after physical activity. By using it before exercise, MIRoTEC helps warm up muscles, reducing the risk of "cold-start" strains.

Furthermore, MIRoTEC proves especially helpful in cold conditions, providing warmth and support. It can be used in conjunction with massage and other natural remedies to enhance their effectiveness. Massaging the injured area between MIRoTEC applications can offer additional benefits.


Overall, MIRoTEC is a versatile and user-friendly physiotherapy aid that can help reduce health problems and expedite recovery time for various musculoskeletal conditions and injuries.



The thermography photos below demonstrate the rise in temperature from use of a Mirotec Back Wrap



Before



After 21 Minutes



Email: sales@gillgroupinternational.com

Mobile: +61 (0)417 388 008

www.gillgroupinternational.com

MIRoTEC offers a range of wraps and inserts designed for specific body parts and purposes. Here are the available options:



Neck Wrap:

Designed to provide support and promote circulation in the neck area.



Elbow Wrap:

Intended for the treatment and support of the elbow joint.



Uni-wrap:

A versatile wrap that can be used for the arm, wrist, ankle, or foot.



Mitten:

Designed to cover and support the hand.



Vest:

Provides support and heat reflection for the upper body.



Shoulder/Arm Wrap:

Targets the shoulder and arm area for support and improved circulation.



Abdominal Wrap:

Designed to aid in the treatment of cramping and abdominal pain.



Lower Back Wrap:

Provides support and heat reflection for the lower back.



Leg Wrap:

Targets the leg area for support and improved circulation.



Knee Wrap:

Designed specifically for the knee joint.



Boot Wrap:

Intended for the treatment and support of the foot and ankle.

Additionally, MIRoTEC offers inserts that can be used with bedding, chairs, and massage equipment to provide enhanced therapeutic effects.

It is important to take certain precautions when using MIRoTEC products:

If excessive sweating occurs, it is recommended to use MIRoTEC for short periods of time, as the porous construction allows water vapor to escape. In cases of bruising, it is advised to apply ice for 36 hours before using MIRoTEC wraps.

It is not recommended to use any other heating agents or liniments in conjunction with MIRoTEC products.